


- SA Shaper Amplifier
- PF Thyatron Pulse Former
- PA Pulse Amplifier
- T Trigger Tube
- D Delay
- DC A DC Amplifier


# SINYPDET

Note: Leads labeled with Command Symbols are connected to TTS or CTS output as the particular command.  
Nos. in paren. identify misc. F.F.'s


## Start Clock 4100 Start Clock Unit




## WAIT-RESUME 4200 WAIT-RESUME UNIT




## CLOCK 4600 Clock generator




## STOP CLOCK PULSED RELAY 4400



## RELEASE LOAD-RUN SELECTION



## Stop after Command exec.



## STOP RELAY CONTROL 4500

